

Highend Data Acquisition Systems


**ZAHNER**

M E S S Y S T E M E

*Made in Germany*

# Scientific Instrumentation

for Photons & Electrochemistry



# ZAHNER ZENNIUM / IM6

## electrochemical workstations

### General

Our R&D team managed to create instruments with outstanding features, state-of-the-art hardware and an advanced software. The result is a milestone in scientific instrumentation. ZENNIUM and IM6 were developed using our 30 years of experience in producing high-precision electrochemical workstations of the high-end class. They provide a frequency range up to 4 MHz (ZENNIUM) / 8 MHz (IM6), an output current up to  $\pm 2.5$  A (ZENNIUM) /  $\pm 3.0$  A (IM6) and fast signal processing. Special measurement techniques guarantee an ultra high accuracy and a minimal interference with the test object.

ZENNIUM/IM6 come bundled with the outstanding THALES-Z (Zennium release) software package which offers all standard methods and more at a mouse click. This is why the ZENNIUM/IM6 can easily be adapted to very different measurement requirements. Furthermore, with the manifold options available, the ZENNIUM/IM6 is able to grow with its tasks. It is best suited for investigations on fuel cells, batteries and solar cells as well as on membranes and sensors or on coatings and laminates, to name only a few.

ZAHNER-elektrik is known to provide competent service all around the world. Our experienced specialists help you to plan, set up and analyze your experiments in electrochemistry, physics, material science and electronics.

### Hardware

#### General


The hardware of the ZENNIUM/IM6 provides

- ultra low-noise potentiostat
- wide frequency range dual DDS FRA
- high CMRR precision U/I-amplifiers
- PulSAR® state-of-the-art differential 18 bit ADCs for AC
- connectors optimized for High Z & Low Z
- 4/9 extension card slots
- 410 MIPS (Dhrystone 2.1) V4e ColdFire® signal processor
- floating USB 2.0 interface

#### Accuracy

The highest priorities for the development and production of the ZENNIUM instruments are accuracy and reliability. The accuracy map of the ZENNIUM clearly shows the high quality of the hardware. These specifications are proofed by measurements for practical use, based on the high-end components we use. They are not only calculated by the theoretical specifications of some components.

PulSAR® is a registered trademark of Analog Devices, Inc.  
ColdFire® is a registered trademark of Motorola, Inc.



At 10mV amplitude.

Option	Function	ext	Int	EPC42 needed
PAD4	4 channel parallel AD converter		X	
TEMP/U	2 inputs for thermocouples + 2 voltage inputs		X	
DA4	4 analog outputs		X	
RMux	Relay multiplexer for the internal potentiostat		X	
PwrMux	Power multiplexer for the PP series potentiostats	X	X	
TR8M	Transient recorder up to 40 MHz		X	
HiZ probe	High impedance probe set	X		
fF probe	Low capacitance probe set	X		
LoZ	Cable set for low impedances	X		
EPC42	Control module for up to 4 external potentiostats		X	
XPot	External standard potentiostat	X		X
PP series	External power potentiostats	X		X
EL series	External high current one quadrant potentiostats	X		X
NProbe	Probe set for measuring electrochemical noise	X		X
COLT	Set-up for coating and laminate testing	X		
CImps	Set-up for photo electrochemical applications	X		X
EiChem Cells	KMZ and AMZ type cells for various applications			



# ZAHNER ZENNIUM / IM6

## electrochemical workstations

### External Potentiostats and Electronic Loads



... enhance the  
output up to  
 $\pm 120$  Volts  
 $\pm 40$  Ampere  
+ 600 Ampere (load)  
50 000 Watt  
64 Channels

### NetVI Interface



... connect the  
Zennium/IM6 easily  
with external data  
acquisition & control  
devices over TCP/IP as  
virtual net instruments

### Up to 9 Extension Card Slots

... expand your Zennium/IM6 optionally  
with plug & play cards for fast data acquisition  
and control output channels according  
to your special requests.

ethernet

### The Modular System

... the only  
electrochemical workstation  
based on a universal modular  
data acquisition system



### Probes

... prepare the Zennium  
for special fields  
of application like  
high impedances,  
electrochemical noise,  
very low capacitance,  
etc. and interface with  
3rd-party devices  
(electronic loads,  
potentiostats, ...)



# ZAHNER ZENNIUM / IM6

## THALES Z software package

### Software

The powerful Thales software package is part of the standard equipment of each IM6 and ZENNIUM system. It provides a multitude of measuring and analysis methods.

The Thales software provides unique features for the acquisition and analysis. The most prominent examples are SCRIPT and SIM. SCRIPT offers you to program user defined combinations of all types of electrochemical measurements, mathematical analysis, algorithms, documentation and data export to a reproducible, fully automatic process. With the outstanding features of SIM you are able to create equivalent circuits and fit the measurement data to these models. The ZHIT feature of SIM helps you to validate your impedance spectra.

For evaluation only it is also possible to run Thales software package on a ThalesBox. The ThalesBox provides a single user licence of Thales. So you can process your recorded data on a simple PC or Laptop while your tests still run on the IM6 or ZENNIUM ...

- System requirements: IBM compatible PC  
Microsoft® Windows® XP / Vista / 7 / 8  
(32-bit / 64-bit)  
USB 2.0

#### General Fields of Application

- low impedance applications  
(fuel cells, batteries, super-caps ...)
- high impedance applications  
(coatings, laminates, membranes, sensors, corrosion ...)
- photoelectrochemical applications  
(silicon, dye-sensitized and organic solar cells,  
organic LED, semiconducting films ...)

#### Look & Feel

- electrochemical methods guide
- free scaleable application window
- Windows® tooltips
- ThalesViewer (Windows® Explorer extensions)
- fast online context sensitive help system
- icon navigation + direct access pull down menus
- child windows for configurable online measurement display

#### Special Functionalities

- fast multiple filter for impedance and photovoltage /current spectra
- special modeling support for solar cells
- joint multiple transfer function fitting
- SCRIPT procedures for user-defined measurement, analysis- and documentation tasks
- multi-channel measuring data acquisition and control in parallel to the electrochemical experiments

#### Connectivity

- remote control via LabVIEW® VI
- integrate third-party signal acquisition devices over TCP/IP as NetVI
- ASCII data logging via online display

Windows®  
tooltips on  
icon touch

Context sensitive  
help system  
on alternate  
mouse button

Fast acquisition  
channels over  
plugin cards

Virtual net  
instruments  
over TCP/IP




# ZAHNER ZENNIUM / IM6

## THALES Z software package



Export of ASCII text tables, vector- and bitmap graphics on all display panels



Direct access pull down menus

Universal Screenshot Function

Remote control with LabVIEW® over TCP/IP anywhere

Output control channels

Graphic and legend preview with Windows® Explorer®

### Supported Method

### Option\*

Testsampling Online Display and Logging	-
Parallel Testsampling Online Display and Logging	PAD4
Electrochemical Impedance Spectroscopy	-
Parallel Electrochemical Impedance Spectroscopy	PAD4
EIS Series Measurements vs. Parameters	-
Impedance Network Analysis, Simulation & Fitting	-
Impedance/Parameter, Capacity/Voltage Curves	-
Parallel Impedance/Parameter, Capacity/Voltage Curves	PAD4
AC Voltammetry	-
Cyclic Voltammetry	-
Linear Sweep Voltammetry	-
Triangular Wave Voltammetry	-
Cyclic Triangular Wave Voltammetry	-
Chronoamperometry With Multiple Triangular Potential Sweeps	-
Stationary Electrode Voltammetry	-
Chronoamperometry With Linear Potential Sweep	-
Tafelscan	-
Steady State Current/Voltage Curves	-
Multi-Cell Current/Voltage Curves	XPOT, PP series
Corrosion / Polarisation Measurement	-
Potentiometry	-
Chronopotentiometry	-
Chronoamperometry	-
Pulse Plating	-
Universal Current/Potential/Time Curves	-
Time Domain Controlled Measurements	-
Electrochemical Noise	NProbe
CorrEISNoise	NProbe
Tast Polarography	-
Differential Pulse Polarography	-
Differential Pulse Voltammetry	-
Stripping Voltammetry	-
Standard Addition Measurement & Analysis	-
Universal Measurement Data Acquisition & Control	NET*
Universal Frequency Response Analysis	-
AC-DC-AC Tests	COLT
Layer Quality Test / Bi-Layer Test	COLT
Fast Pulse & Transient Recording	TR8M
High Current Interrupt Measurements	HCI
Relaxation Voltammetry	NProbe
Solar Cell Fill-Factor, Efficiency, Maximum Power, OCP, ISC	CIMPS
Controlled Intensity Modulated Photocurrent Spectroscopy	CIMPS
Controlled Intensity Modulated Photovoltage Spectroscopy	CIMPS
Charge Extraction After N. W. Duffy, L. M. Peter et. al.	CIMPS
Light Transient Measurements. Chopped Light Voltammetry	CIMPS
DC vs. Intensity Transfer Functions, Time Domain Measurements	CIMPS
Lightemission Voltage/Current Characteristic (OPV)	CIMPS +
Spectral Resolved OPV, Absolute PL Quantum Yield	CIMPS +
Spectral Absorption Voltage/Current Characteristic	CIMPS +
Photocurrent Spectroscopy	CIMPS +
Quantum Efficiency, Incident Photon to Charge Carrier Efficiency	CIMPS +
Photocurrent/- voltage Response on Fast Light Transients	CIMPS +
Dynamic and Static Transmittance / Reflectance	CIMPS +
Synchronous Multi Spectral DTR with Parallel Impedance	CIMPS +
Cross Transfer Function Measurements	-
Membrane Penetration Transfer Function Measurements	-
Programmable Procedures	-

\* Required hardware options

# ZAHNER ZENNIUM/ IM6

## electronic load EL1000



**EL1000**

### Electronic Load EL1000

Electronic loads are indispensable tools in several fields of electrochemistry, for example in the research of batteries and fuel cells. The EL 1000 was designed to investigate single cells as well as complete stacks, either as stand-alone device under PC control or in combination with a workstation IM6 or Zennium for instance for impedance measurements. Adding an additional external electronic load, the power can be raised up – adding the PAD4 to the controlling workstation, individual segments of the stack can be investigated synchronously in parallel.

#### Specifications


Frequency range	10 $\mu$ Hz - 100 kHz
Current range	200 A / 600 A (with 3 <sup>rd</sup> party electronic load)
Voltage range	$\pm 4$ V / $\pm 100$ V
Maximum power	1.000 W (stand-alone) scaleable with 3 <sup>rd</sup> party electronic load
Dimensions	470 x 160 x 446 mm

**Optional:**  
electrochemical workstation  
IM6/Zennium  
for dynamic measurements

- impedance spectra (EIS)
- series EIS



**Optional:**  
additional electronic load  
for high power (< 600 A)





# ZAHNER ZENNIUM / IM6

## parallel impedance add-on

### PAD4: 4 Channel Synchronous Impedance Converter

#### True Parallel Synchronous Impedance

**Save time - measure up to 17 stack-cells in one run - no time mismatch between impedance spectra - record anodic, cathodic & total impedance simultaneously - measure additional transfer function signals...**

The Zahner PAD4 is a 4-channel add-on card for Zahner Electrochemical Workstations. It introduces four additional parallel sampled signal inputs for cell voltage and impedance in fuel cell stacks and battery packs, with a common current. The Zahner Zennium supports up to two PAD4-cards for up to nine parallel channels, while the IM6 can control up to four cards for a maximum count of 17 parallel channels.


PAD4-cards are plug 'n' play – they are detected automatically on start-up. The PAD4 may be combined with the basic ECW or with the ECW controlling a slave potentiostat or an electronic load, finally providing tests on stacks of up to 100V / 600A / 50 KW.

#### Additional Methods


- parallel electrochemical impedance spectroscopy
- parallel impedance/parameter capacitance/voltage curves
- parallel testsampling online display and logging

#### Specifications

Channels / card	4 individually addressable
Impedance measurement:	
Frequency range	10 $\mu$ Hz to 250 kHz
DC-potential measurement:	
Voltage range	$\pm 4$ V
Common mode range	$\pm 100$ V
A/D converter resolution	18 bit



PAD4 Nyquist plot of a five cell SOFC stack



Up to four  
PAD4 add-on cards  
allow the parallel measurement  
of up to 17 cells of a stack

# ZAHNER ZENNIUM / IM6

## low capacitance probe



### fF-Probe: femto-Farad Probe

#### Low Capacitance Measurement Probe Set

The femto-Farad Probe works as a front-end to the IM6/Zennium potentiostat. Apart from its limited current capability, all basic functionalities of the Thales software are supported. In particular impedance spectroscopy can be applied. Due to the fact, that the primary measurement magnitude is the complex impedance, besides the sample capacity, resistive and DC contributions can be determined as well.


#### Additional Methods

- low capacitance electrochemical impedance spectroscopy
- low capacitance impedance/parameter

#### Specifications

Frequency range	10 $\mu$ Hz to 1 MHz
Current auto ranging, defeatable current ranges	0nA - $\pm 40$ nA $\pm 40$ nA - $\pm 400$ nA $\pm 400$ nA - $\pm 4$ $\mu$ A $\pm 4$ $\mu$ A - $\pm 40$ $\mu$ A
Voltage range	$\pm 4$ V
Resolution of any range	18 bit
Capacity offset	$\pm 1$ fF <sup>*)</sup>
Capacity resolution	$\pm 0,1$ fF <sup>*)</sup>
Capacity accuracy	$\pm 0,25\%$ of reading $\pm 2$ fF <sup>*)</sup>


<sup>\*)</sup> current range  $\pm 40$  nA,  
AC amplitude  $\geq 100$  mV,  
zero DC current




The trans-impedance principle  
for the determination of small  
capacities

SO "hot" test signal output  
VGI "virtual ground" signal input  
Z device under test  
V, A voltage and current measurement  
C<sub>si</sub> and C<sub>so</sub> parasitic stray capacitance at the input and output terminal

Microelectromechanical systems  
(MEMS) and Sensors



Determination of the  
coupling capacity between two  
adjacent pads of a printed circuit board



Nanoelectrodes



# ZAHNER ZENNIUM / IM6

## PP-Series power potentiostats

### PP2x1 Power Potentiostats

The PP-Series potentiostats are designed to apply and sink high currents up to  $\pm 40\text{A}$  at a total power dissipation of up to 200W.

The PP-Series potentiostats are controlled by an EPC42, a plug-in module for the Electrochemical Workstations IM6 and Zennium.

Up to four EPC42 cards can be installed in an IM6 or an Zennium. In total, up to 16 PPs can be controlled by one IM6/Zennium. Each potentiostat will hold the control

parameters from one access to the next one, so that no potential or current disturbances can occur while scanning the potentiostats. If series measurements are performed with more than one PP-Series potentiostat, spectra are taken from all modules in a definable order.

The PP-Series is embedded completely in the IM6/Zennium environment. Thus, all acquisition and analysis techniques that run on the IM6/Zennium can be applied with the power potentiostats as well. The installation of one or more PP-Series potentiostat will upgrade your IM6/Zennium to an even more versatile, high-current electrochemical workstation.

The PP-Series potentiostats can also be controlled by a Windows®-PC. In this case they provide methods, summarized in the table below. They also work as a LabVIEW™ Virtual Instrument under the LabVIEW™ software. To implement the PP-Series potentiostats into existing test environment, a supporting DLL is available on demand.

You can control the PP-Series potentiostats in a mixed mode with an IM6/Zennium and a PC in parallel. Both devices can be connected and disconnected during operation.



#### Supported Methods with IM6/Zennium

- impedance spectroscopy
- simulation & model fitting
- cyclic voltammetry
- polarisation curves
- multicell multitasking voltammetry
- arbitrary current/potential/time measurements
- capacity/potential measurements
- automatic series measurements

#### Supported Methods with PC

- test sampling
- U vs. time, I vs. time
- current potential curves (U/I)
- cyclic voltammetry
- charging/discharging, battery cycling
- LabView® virtual instrument
- DLL support available

#### Specifications

Model name	PP201	PP211	PP241
Operating modes	pot/gal/oc	pot/gal/oc	pot/gal/oc
Potential range	$\pm 10\text{ V}$	$\pm 20\text{ V}$	$\pm 5\text{ V}$
Potential accuracy	$\pm 0.1\% / \pm 1\text{ mV}$	$\pm 0.1\% / \pm 2\text{ mV}$	$\pm 0.1\% / \pm 1\text{ mV}$
Current range	$0\text{ A} \dots \pm 20\text{ A}$	$0\text{ A} \dots \pm 10\text{ A}$	$0\text{ A} \dots \pm 40\text{ A}$
Current accuracy	$\pm 0.25\% / \pm 1\text{ mA}$	$\pm 0.25\% / \pm 1\text{ mA}$	$\pm 0.25\% / \pm 1\text{ mA}$
Output power	200 W	200 W	200 W
Frequency range	10 $\mu\text{Hz}$ - 200 KHz	10 $\mu\text{Hz}$ - 200 KHz	10 $\mu\text{Hz}$ - 200 KHz
Impedance range	1 $\mu\Omega$ - 1 K $\Omega$	1 $\mu\Omega$ - 1 K $\Omega$	1 $\mu\Omega$ - 1 K $\Omega$
Ambient temperature	0 °C ... 30 °C	0 °C ... 30 °C	0 °C ... 30 °C
System requirements	IM6/Zennium+EPC42 or PC	IM6/Zennium+EPC42 or PC	IM6/Zennium+EPC42 or PC

# ZAHNER ZENNIUM / IM6

## specifications

General	Zennium	IM6
Overall Bandwidth	DC - 5 MHz	DC - 10 MHz
ADC Resolution	18 bit	
Harmonic Reject	> 60 dB @ 1/2 full scale	
Potentiostat Modes	Potentiostatic, galvanostatic, pseudo-galvanostatic, rest potential, off	
Cell Connection	2-, 3-, 4-terminal Kelvin	
Chassis	ground	
Extension Slots	4	9 (incl. 1x EPC42)
PC interface	USB 2.0	
Dimensions	364 x 160 x 376 mm	470 x 160 x 376 mm
Weight	12 kg	15 kg
Accessories	U-buffer, 2 cell cable set, USB-cable, power cord, manual	+Thalesbox, +EPC42
Power supply	230/115 V, 50/60 Hz	
Ambient temperature	+10° C to +30° C	
Ambient Humidity	< 60% without derating	

### Frequency Generator & Analyzer

Frequency Range	10 $\mu$ Hz to 4 MHz	10 $\mu$ Hz to 8 MHz
Accuracy	< 0.0025%	
Resolution	0.0025%, 10.000 steps/decade	

### Output Potentiostatic

Controlled Voltage	Pot $\pm 4$ V U-buffer $\pm 10$ V	
Resolution	Pot 125 $\mu$ V U-buffer 320 $\mu$ V	
Accuracy	Pot $\pm 250$ $\mu$ V U-buffer $\pm 1$ mV	
Temperature Stability	better 20 $\mu$ V/°C	
Compliance Voltage	Pot $\pm 14$ V CVB120 $\pm 120$ V	
AC-Amplitude	1 mV to 1 V	
Bandwidth	4 MHz @ 33 $\Omega$ load	8 MHz @ 33 $\Omega$ load
IR Compensation	Method Auto AC Impedance Technique Range 0 to 10 M $\Omega$ Resolution 0.012%	
Small Signal Rise Time	250 ns to 200 $\mu$ s in 5 steps, automatic selection	
Slew Rate	15 MV/s	
Phase Shift	10° @ 250 kHz	

### Output Galvanostatic

Controlled current	$\pm 2.5$ A	$\pm 3.0$ A
Current Range	Pot $\pm 100$ nA to $\pm 2.5$ A in 10 steps HiZ $\pm 1$ nA to $\pm 0.5$ A in 12 steps	$\pm 100$ nA to $\pm 3.0$ A in 10 steps
Min. Resolution	0.025%	
Accuracy	0.1% @ > 2 $\mu$ A to 100 mA 1% @ < 2 $\mu$ A or > 100 mA	

### Input

Potential Ranges	Pot $\pm 1, \pm 2, \pm 4$ V U-buffer $\pm 4, \pm 10$ V	
Potential Resolution	256.000 steps per range	
Offset Voltage	< 100 $\mu$ V	
Offset Temperature Stability	< 20 $\mu$ V/°C	
Current Range	Pot $\pm 100$ nA to $\pm 2.5$ A in 30 steps, automatic range selection HiZ $\pm 1$ nA to $\pm 0.5$ A in 32 steps, automatic range selection	$\pm 100$ nA to $\pm 3.0$ A
Accuracy	0.05% @ > 2 $\mu$ A to 100 mA 0.5% @ < 2 $\mu$ A or > 100 mA	
Input Bias Current	Pot 10 pA HiZ 12 fA	
Current Resolution	Pot 2.5 pA HiZ 25 fA	
Input Impedance	Pot 1 T $\Omega$ // $\pm 5$ pF (typical), 50 G $\Omega$ // $\pm 5$ pF (min.) HiZ 10 T $\Omega$ // $\pm 1$ pF (typical), 1 T $\Omega$ // $\pm 1$ pF (min.)	
Impedance Range	Pot 1 m $\Omega$ to 1 G $\Omega$ / 2% 100 m $\Omega$ to 10 M $\Omega$ / 0.2% HiZ 100 m $\Omega$ to 100 G $\Omega$ / 3% Gal 30 $\mu$ $\Omega$ to 1 G $\Omega$ / 2%	
Common Mode Rejection	> 86 dB @ 10 $\mu$ Hz to 100 kHz > 66 dB @ 100 kHz to 4 MHz	> 66 dB @ 100 kHz to 8 MHz
Input Channel Phase-Tracking accuracy	$\pm 0.1^\circ$ @ 10 $\mu$ Hz to 100 kHz $\pm 0.25^\circ$ @ 100 kHz to 4 MHz	$\pm 0.25^\circ$ @ 100 kHz to 8 MHz
Equiv. Effective Input Noise	20 $\mu$ V rms / 2 pA rms @ 1 mHz to 10 Hz	